

Implant compatibility chart



Manufacturer	Implant	Titanium GoldHue	Zirconia, including ATLANTIS™ crown abutment	Screwdriver
BioHorizons	INTERNAL			
	Internal 3.5, 4.5, 5.7	x	x	1.25 mm Hex screwdriver
	Tapered Internal 3.5, 4.5, 5.7	x	x	
	EXTERNAL			
	External Implant 3.5	x		0.050" Hex screwdriver
External Implant 4.0, 5.0	x	x		
Biomet 3i	INTERNAL			
	Certain™ 3.25, 4.0, 4/3, 5.0, 5/4, 6.0	x	x	0.048" Hex screwdriver
	Certain™ Prevail™ 3/4/3 (3.4), 4/5/4, 5/6/5, 4/3, 5/4, 6/5	x	x	
	Certain™ XP 4/5, 5/6	x	x	
	EXTERNAL			
	MicroMiniplant™ 3.25	x		0.048" Hex screwdriver
	Miniplant™ 3.25	x	x	
	Osseotite™ 3.75, 4.0, 5.0, 6.0	x	x	
Osseotite™ XP 4/5, 5/6	x	x		
Camlog	INTERNAL			
	Screw-line Implant 3.3	x		1.28 mm Hex screwdriver
	Screw-line/Root-line Implant 3.8, 4.3, 5.0, 6.0	x	x	
DENTSPLY Implants <small>Note! OsseoSpeed™ implants are part of the ASTRA TECH Implant System™</small>	INTERNAL			
	OsseoSpeed™/OsseoSpeed™ TX 3.0 – Yellow	x		0.050" Hex screwdriver
	OsseoSpeed™/OsseoSpeed™ TX 3.5, 4.0 – Aqua	x	x	
	OsseoSpeed™/OsseoSpeed™ TX 4.5, 5.0 – Lilac	x	x	
	OsseoSpeed™ TX Profile 4.5, 5.0	x	x	
	OsseoSpeed™ 4.0 – Blue (Previous)	x		1.22 mm Hex screwdriver
	XiVE® S 3.0	x		
	XiVE® S 3.4, 3.8, 4.5, 5.5	x	x	
Frialit™ 3.4, 3.8, 4.5, 5.5	x	x		
Keystone Dental	INTERNAL			
	PrimaConnex™ SD 3.3, 3.5, Genesis™ 3.8	x	x	Keystone Quad Driver
	PrimaConnex™ RD 4.0, 4.1, Genesis™ 4.5	x	x	
	PrimaConnex™ WD 5.0, Genesis™ 5.5, 6.5	x	x	
	EXTERNAL			
	Restore™ SD 3.3, 3.4	x		0.048" Hex screwdriver
	Restore™ RD 3.75, 4.0	x	x	
Restore™ WD 5.0, 5.5, 6.0	x			
Nobel Biocare	INTERNAL			
	NobelActive™ NP 3.5	x	x	Nobel Biocare Unigrip
	NobelActive™ RP 4.3, 5.0	x	x	
	NobelReplace™ NP 3.5	x	x	
	NobelReplace™ RP 4.3	x	x	
	NobelReplace™ WP 5.0	x	x	
	NobelReplace™ 6.0	x	x	
	EXTERNAL			
	Brånemark System™ NP 3.3	x		Nobel Biocare Unigrip
	Brånemark System™ RP 3.75, 4.0	x	x	
	Brånemark System™ WP 5.0, 5.5	x	x	
	SteriOss Replace™ 3.1, 3.5, 4.3, 5.0, 6.0	x		0.050" Hex screwdriver
SteriOss HL™ 3.25, 3.8, 4.5, 5.0, 6.0	x			
Straumann	INTERNAL			
	Bone Level 3.3 NC (Narrow CrossFit™)	x		SCS screwdriver
	Bone Level 4.1, 4.8 RC (Regular CrossFit™)	x	x	
	Standard/Standard Plus 4.8 RN (Regular Neck)	x		
	Standard/Standard Plus 4.8 WN (Wide Neck)	x		
	EXTERNAL			
Standard Plus 3.3 NN (Narrow Neck)	x		SCS screwdriver	
Zimmer Dental	INTERNAL			
	Tapered Screw-Vent™ 3.5, 4.5, 5.7	x	x	1.25 mm Hex screwdriver
	Screw Vent™ 3.3, 3.7, 4.5	x	x	
	SwissPlus™ 4.8	x		

All trademarks are the property of their respective owners.